




Golden Star Surgical Industries
MANUFACTURER & EXPORTER

Manual For use of Castraion

Prevention:

Hygiene and prophylactic use of vaccines can help reduce complications from castration. Important considerations include:

- Use disinfectants such as iodine to clean the scrotum prior to castration
- Prophylactic antimicrobials can help reduce infection at surgical or banding site
- Tetanus toxoid and vaccination can reduce the incidence of tetanus infection due to open wounds from castration
- Provide a clean, dry environment for calves after castration to avoid infections
- Monitor calves closely for two weeks after castration
- Look for swelling, signs of infection, tetanus and abnormal gait
- Treat wounds as needed
- Time castration properly to avoid fly season. Open wounds from castration during fly season can lead to flystrike and infection.



By www.goldenstarsurgical.com for video results please visit <http://www.youtube.com/watch?v=SzY-xZDt9P0>



Golden Star Surgical Industries

MANUFACTURER & EXPORTER

Product Picture cum catalogue for Burdizzo Castrator 19" Imported

Material- Carbon steel, min. 0.4% carbon, heat treated to 30-40 HRC, length-460-480 MM, Upper jaw width-75mm+/- 2mm, lower jaw- stopcord 5-8mm ht on both sides, Gap on either side between upper jaw and lower jaw should range between 1-1.5 mm.

Along With Carry Bag.





Golden Star Surgical Industries
MANUFACTURER & EXPORTER

Product Picture cum catalogue for Burdizzo Castrator 9" Imported

Material- Carbon steel, min. 0.4% carbon, heat treated to 30-40 HRC,length-220-230 mm,Upper jaw width-42mm+/- 2mm, lower jaw- stop cord 5-8mm ht on both sides,Gap on either side between upper jaw and lowerjaw should range between 1-1.5 mm.

Along with Black Carry Bag



Visit us on www.goldenstarsurgical.com A Name you can trust Upon Since 1985